

EEG Acquisition



Preprocess


1-30Hz filter

baseline
correction

50Hz notch
filter

Downsample

Offline Training




Best Model Determination

Window
selection

Model
Optimization

EEG Acquisition

Real-time Data Stream



Online EEG Acquisition



Real-time Decoding Module

&

Data
Communication

Synchronization

Send Decoding
Results

Preprocess

Model
Depoly

Model
Decoding

Receive
Command

Real-time Interaction Module

Data
Communication

Receive
Results

Send Receive
Results

User
Interface

Unity



Feedback